

Heating Bath



Heating baths are used in a variety of laboratory environments to provide a precisely controlled temperature that remains constant throughout testing. We offer several sizes of the Heating Bath, please see following page for full range.

Key Features

- Stainless steel interior
- Digital display, Microprocessor control with timing function.
- Audible and visible alarm indicates if over temperature
- Power Supply : 220V 50HZ
- Temperature Stability : ± 0.5 °C
- Temperature Track Alarm : +2 °C
- Timing Range : 1~9999min

Ascott Analytical Equipment Limited

6-8 Gerard, Lichfield Road Industrial Estate, Tamworth, Staffordshire, B79 7UW, Great Britain

T +44 (0) 1827 318040 F +44 (0) 1827 318049 E info@ascottshop.com W www.ascottshop.com

Available Models

Product Code	Temperature Range	Power	Cubage	Interior Dimensions (mm)W×D×H	Overall Dimensions (mm)W×D×H
T01-701131	RT+5~99 °C	500W	11L	420×180×150	570×270×260
T01-701132	RT+5~99 °C	1000W	34L	600×300×190	750×400×300
T01-701133	RT+5~99 °C	400W	12L	320×240×160	460×280×190
T01-701134	RT+5~99 °C	600W	22L	450×300×160	610×340×190
T01-701135	RT+5~99 °C	1000W	30L	600×300×160	760×340×190



Ascott Analytical Equipment Limited

6-8 Gerard, Lichfield Road Industrial Estate, Tamworth, Staffordshire, B79 7UW, Great Britain

T +44 (0) 1827 318040 F +44 (0) 1827 318049 E info@ascottshop.com W www.ascottshop.com