

Low Temperature Thermostatic Bath



Vertical



Horizontal

Low Temperature Thermostatic Baths are widely used in the fields of biological engineering medicine food chemical industry chemical analysis petroleum, etc., as they provide precise controlled temperatures. They are very responsive and can heat and cool quickly and maintain constant temperatures with a high degree of accuracy.

See following page for the full range of models that we offer

Ascott Analytical Equipment Limited

6-8 Gerard, Lichfield Road Industrial Estate, Tamworth, Staffordshire, B79 7UW, Great Britain

T +44 (0) 1827 318040 F +44 (0) 1827 318049 E info@ascottshop.com W www.ascottshop.com

Available Models of Low Temperature Thermostatic Bath

Product Code	Range of Temp (°C)	Fluctuating degree of Temp (°C)	Volume of inner container (mm)	Volume of inner	Size of opening (mm)	Depth (mm)	Pump flux (L/min)
T01-701142	0~100	±0.05	260×170×120	6 L	150×150	120	6
T01-701143	-5~100	±0.05	250×200×150	7.5 L	180×150	150	6
T01-701144	-5~100	±0.05	300×250×200	15 L	235×180	200	6
T01-701145	-5~100	±0.05	400×325×230	30 L	310×280	230	13
T01-701146	-10~100	±0.05	250×200×150	7.5 L	180×150	150	6
T01-701147	-20~100	±0.05	250×200×150	7.5 L	180×150	150	6
T01-701148	-20~100	±0.05	300×250×200	15 L	235×180	200	6
T01-701149	-20~100	±0.05	400×325×230	30 L	310×280	230	13
T01-701150	-30~100	±0.1	250×200×120	6 L	180×150	120	6
T01-701151	-30~100	±0.1	250×200×200	10 L	180×150	200	6
T01-701152	-30~100	±0.1	300×250×200	15 L	235×180	200	6
T01-701153	-30~100	±0.1	400×325×230	30 L	310×280	230	13
T01-701155	0~100	±0.05	290×180×150	8 L	155×155	150	8
T01-701156	-10~100	±0.05	290×180×150	8 L	155×155	150	6
T01-701157	-21~100	±0.05	280×250×130	9 L	235×180	130	6



Ascott Analytical Equipment Limited

6-8 Gerard, Lichfield Road Industrial Estate, Tamworth, Staffordshire, B79 7UW, Great Britain

T +44 (0) 1827 318040 F +44 (0) 1827 318049 E info@ascottshop.com W www.ascottshop.com