

Viscosity Calibration Oils

| Ordering Information | Kinetic Viscosity cP | Kinematic Viscosity mm ² /s | Suitable for |
|----------------------|-------------------------|---|--|
| T01-700336 | 5 | 5.4 | ISO Flow Cup No.3 |
| T01-700337 | 50 | 52.1 | Zahn Cup No.1 Ford Cup No.2 |
| T01-700338 | 100 | 104.2 | Zahn Cup No.2 ISO Flow Cup No.4 Ford Cup No.3 Iwata Cup No.2 NK-2 |
| T01-700339 | 200 | 208.3 | ISO Flow Cup No.5 DIN Flow Cup No.4 Ford Cup No.4 |
| T01-700340 | 500/67.2 KU | 520.8 | Zahn Cup No.3 and No.4 ISO Flow Cup No.6 |
| T01-700341 | 1,000/86.3 KU | 1,052.3 | Zahn Cup No.5 Krebs Stormer Viscometer |
| T01-700342 | 5,000 | 5,170.6 | Rotary Viscometer |
| T01-700343 | 10,000 | 10,416.7 | Rotary Viscometer |
| T01-700344 | 30,000 | 31,039.8 | Rotary Viscometer |
| T01-700345 | 60,000 | 62,176.2 | Rotary Viscometer |
| T01-700346 | 100,000 | 103,092.8 | Rotary Viscometer whose range is over 100,000 cp |
| T01-700347 | 67.2 KU | ---- | Krebs Stormer Viscometer |
| T01-700348 | 86.3 KU | ---- | Krebs Stormer Viscometer |
| T01-700349 | 106.6 KU | ---- | Krebs Stormer Viscometer |
| T01-700350 | 121.5 KU | ---- | Krebs Stormer Viscometer |

Ascott Analytical Equipment Limited

6-8 Gerard, Lichfield Road Industrial Estate, Tamworth, Staffordshire, B79 7UW, Great Britain

T +44 (0) 1827 318040 F +44 (0) 1827 318049 E sales@ascott-analytical.com W www.ascottshop.com